

Palium

Automatic reusable devices for soft tissues biopsy

MDL



Palium (PL) is an automatic, reusable device equipped with a sequential advancement system of cannula and stylet.

Designed to get a high quality of the core specimen.

Ultra-sharp tips and polished surfaces for a better entry and exit.

Powerful **dual spring system** for consistent and excellent core samples.

Rapid-fire action.

Small size and lightweight. One hand activation to provide an easy access to the tissue sample.

Dual penetration depths for flexibility and convenience **(22mm or 15mm)**.

The samples can be taken with just one shot of the button.

Can be easily **wielded with just one hand**. May be sterilized with conventional methods.

FEATURES:

- Side entry for **MGK1 needle** spacer
- Safety to avoid accidental firing
- Trigger site on the rear
- Safety button

Sequential loading system

that requires a slight activation force. Advancement length: 15 - 22mm.



Use exclusively with MDL needle,
Palium needle
(code MGk1)



M.D.L. Srl
Via Tavani 1/a - 23014 - Delebio (So) Italy
Tel +39 -0342 682130 Fax +39 0342 691316
email : uc@mdsrl.com
web : www.mdsrl.com